

Exp 1: Scientific Measurements - CHEM 117-506

Kevin Lei

Lab Partner: Eli Gettleman

9/12/2022

Purpose: To become familiar with the proper procedure of taking scientific measurements and working with lab equipment

Procedure:

- ① Acquire one of each of the following: 125 mL Erlenmeyer flask, 100 mL beaker, and ~~100~~¹⁰ mL graduated cylinder
- ② Pick one of the pieces of glassware and write down which was chosen
- ③ Go to a top-loading balance, tare, and weigh the dry glassware. Record all digits of mass.
- ④ Go to an analytical balance, close the doors, tare, open the doors, then place dry glassware on it. Close doors and record all digits of mass.
- ⑤ If beaker or graduated cylinder: Fill glassware about $\frac{1}{3}$ full with water and record volume. If erlenmeyer flask, using a buret, record initial volume of water in buret, fill flask approximately $\frac{1}{3}$ with water and record final volume of water in buret. (Make sure measurements have one more level of significance than is labeled on glassware)
- ⑥ Weigh the glassware with water on both top-loading and analytical balances and record respective masses
- ⑦ Dump water into the sink weigh the wet glassware on both top-loading and analytical balances and record their respective masses
- ⑧ Repeat steps 2-7 for each remaining piece of glassware
- ⑨ Return all glassware where it was found and wipe down the benchtop with a wet paper towel

Safety: Approved chemical splash goggles must be worn when chemicals and/or glassware are in use anywhere in the lab

Waste Disposal: water is dumped in the sink

Data and Observations

Erlenmeyer flask dry weight

- top-loading: 92.82 g
- analytical: 92.8468 g

Initial volume of water in buret: ~~42.20~~^{0.00} mL

Final volume of water in buret: 50.00 mL

Erlenmeyer flask weight with water

- top-loading: 142.46 g
- analytical: 142.5196 g

Wet Erlenmeyer flask weight after dumping water

- top-loading: 93.63 g
- analytical: 93.6548 g

Beaker dry weight

- top-loading: 67.59 g
- analytical: 67.5817 g

Volume of water in beaker: 40.0 mL

Beaker weight with water

- top-loading: 104.29 g
- analytical: 104.2629 g

Wet beaker weight after dumping water:

- top-loading: 67.79 g
- analytical: 67.8134 g

Graduated cylinder dry weight:

- top-loading: 94.89 g
- analytical: 94.9264 g

Volume of water in graduated cylinder: ~~20.00~~ 20.00 mL

Graduated cylinder weight with water:

- top-loading: 114.01 g
- analytical: 113.9890 g

Graduated cylinder weight after dumping water:

- top-loading: 95.65 g
- analytical: 95.6580 g

W. W. W.